

Ultra Pure Water System

300GPD RO

Lab water system

Optimal system for laboratory and medical use.

- Fully meet the demand water quality.
- Advanced RO system with Deionization system.
- Full automatic control.
- Easy maintenance.



Applications

- Pharmaceutical industry
- Analytical Research Laboratories
- Drug Screening Laboratories
- Clinical Testing Laboratories
- Biotechnology Industry

2-WAY water outlet
&
Dual TDS Monitor



300GPD RO

Water treatment System



600GPD RO



*Optional
Filter System!*

*20inch filter Housing
Bigger Capacity!*



900GPD RO



waco

MADE IN KOREA

the infinite
waco

Commercial Water treatment System



20ton / Day

WACO Corp.
(OLD : Hyundai Wacotec Co., Ltd.)

MADE IN KOREA

Industrial RO Water treatment System

Enough Capacity
5~30 ton



PROCESS FLOW

- Raw Water
- Multi-media iron removal filter
- De-ionizing filter
- Micro filter
- RO Membrane
- Activated carbon block filter
- UV sterilizer
- Smart filter (UF)
- Drinking Water



POE(Point-of-Entry)

To supply clean and safe water to the whole house & restaurant, hospital and etc.

Easy Installation and Operation

Reverse Osmosis System

Portable Bottled Water Factory

Full automatic system managed by one man,
Suitable for small sized water factory, Hotel,
Resort, Military,



- Applied the HEPA filter air purification system for maintaining Clean and hygienic production space.
- Applied the UV filter and UV lamp for sterilization of mineral water and bottle cap to prevent secondary pollutant.
- Available to choose Ozone water system(option) for sterilization of bottle.

R/O water purification system(option)

